**You are not alone!**  
**Join Us!**

A community of women  
for support and advice.

   @IntimateRose

Dr. Amanda Olson, DPT, PRPC

**VELVET ROSE**  
Personal Lubricant

**The Essentials**  
Perfect with the Pelvic Wand


**ENCHANTED ROSE**  
Organic Feminine Balm

**Shop @:**  
**IntimateRose.com**

Manufactured For: Intimate Rose  
North Kansas City, MO 64116  
Made in PRC

 Intimate Rose®

Temperature Therapy  
**Pelvic Wand**

Includes French, German, Italian, Polish, Portuguese,  
Spanish, Swedish, Dutch, and Turkish Translations.

### Indications for Use

The Pelvic Wand is a device developed to manage trigger points in the pelvic floor muscles for the managing of pelvic pain conditions. Such conditions may be the result of injury, trauma, surgery, or chronic muscular conditions such as vaginismus.

### Contraindications

- Pregnant patients should not use the wand unless cleared by their birth provider.
- Patients with active infection of the pelvis or genitals should not use the wand until infection is cleared entirely.

The Temperature Therapy Wand assists in releasing painful trigger points in the pelvic floor using soothing warm or cold temperature and gentle pressure.

To decrease pelvic pain due to tender points, also referred to as trigger points, the keys to success are the consistent use of the Temperature Therapy Pelvic Wand and the practice of pelvic floor muscle relaxation. Your health care provider may have a unique training plan for you outside of the recommendations made here. Always consult with a health care provider before starting a new training plan.

### How to Start

1. For cold therapy, place the wand in the refrigerator or freezer. Prior to use, wash the wand with cold water and mild soap then pat dry.

2. For heat therapy, run the wand under warm water for 1-3 minutes to warm the wand. Finish by gently washing with warm water and mild soap. Do not boil, microwave, or use any other method to warm the wand.
3. Designate a place in your home that is safe, quiet, and comfortable, such as your bed.
4. Decide which end of the wand to use. Either end can be used in the vagina or rectum. The pointed end can be helpful in relieving tender points near the entrance of the vagina, while the longer end can be more helpful in reaching the deeper pelvic muscles.
5. Use a generous amount of water-based lubricant on the first 1 to 2 inches of the desired treatment end of the wand, as well as the opening of the vagina or rectum.
6. Lie on your back with knees bent and feet planted. Some people may prefer to lie on their side instead. To do this, bend your knees and support your top leg with a folded pillow between your knees.
7. Begin by breathing in and allowing your belly to expand, followed by exhaling, allowing your belly to slowly fall. The act of slowly exhaling helps to naturally relax the pelvic floor muscles. Gently bring the wand to the opening of the vagina or rectum and carefully insert it.
8. After insertion into the vagina or rectum, gently press the wand into the internal tissue and move the wand clockwise until encountering a tender point. When you find a tender point, gently compress the end of the wand into the tender point with the same firmness you would use to check a tomato for ripeness. In other words, don't press so hard you squish your tomato.
9. Maintain gentle pressure on the tender point and slowly move your bent knee left and right until you find a position that stops the pain in the pelvic floor muscle. When you find this position, remain there for 1 to 2 minutes to

allow the tender point to fully release. Continue to breathe deeply.  
10. Repeat this process 1 or 2 times per day as needed.

### **General Care and Maintenance**

It is important to clean the wand with water and soap before and after each use. Rinse the wand thoroughly and examine it for any signs of cracks in the medical grade silicone. Small openings near the charging nodes will not cause issues as the wand is water tight. Cracks or openings elsewhere will require a new wand. If you find a crack in your wand, contact us at [Support@IntimateRose.com](mailto:Support@IntimateRose.com).

### **Pain With Penetration**

Pelvic pain due to trigger points is often associated with painful sexual intercourse. Many people have difficulty with penetration and benefit from the use of vaginal dilators, also referred to as vaginal trainers.

Go to [IntimateRose.com/dilators](https://www.intimaterose.com/dilators) to learn more.

### **Warning:**

Non-sterile device is designed for single patient use only. Not manufactured with Natural Rubber Latex.

For use under the recommendation of a health care provider. Do not use silicone-based lubricants. Use Intimate Rose's Velvet Rose lubricant. Do not force the wand into the vagina or rectum. Use slow, gentle movements.

**Indications:** The wand may be used to massage tender points in the pelvic floor via the vagina or rectum.

**Contraindications:** Do not use if you have an acute urinary or genital tract

infection. If heavy bleeding occurs with use of the wand contact your physician.

**Important Note on Hygiene:** Wash hands and wand with mild soap and water prior to use. Rinse thoroughly and pat dry. Apply a generous amount of Intimate Rose Velvet Rose Lubricant to the desired end of the wand and insert slowly into the vagina or rectum.

### **Tender or Trigger Points**

Common symptoms of pelvic pain due to trigger points or tender points include:

- Achy sensations from deep inside the pelvis.
- Burning or zinging sensations during intercourse.
- Burning sensations while urinating without a bacterial infection.
- Sexual pain with no known cause.
- Inability to tolerate vaginal penetration for sex, medical exams or tampon use.
- Pain in the tailbone (coccyx).
- Chronic constipation with an unknown cause.
- Pain in the groin.
- Pain around the bladder.

Trigger points are defined as areas of muscle that are painful to the touch (palpation). They are characterized by the presence of taut bands or “knots” in the muscle and the generation of a referral pattern of pain. Tender points are commonly present in conditions of pelvic pain including vaginismus, vestibulodynia, vulvodynia, coccygeus, proctalgia fugax, endometriosis, post-childbirth, and after pelvic or abdominal surgery.

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